



# Microarrays

## EMD Team Fact Sheet—November 2011

*This fact sheet, developed by the ITRC Environmental Molecular Diagnostics (EMD) Team, is one of 10 designed to provide introductory information about and promote awareness of EMDs. Please review the Introduction to EMDs Fact Sheet along with this one. A glossary is included at the end of this fact sheet.*

### Why are microarrays relevant?

Microarrays offer the ability to simultaneously detect and semiquantitatively assess the relative abundance of thousands of different microbial biomarker genes as a comprehensive evaluation of the microbial community composition and its potential activity within an environmental sample. Microarray analysis offers advantages at sites that require a comprehensive view of the microbial community and where a larger number of biomarker gene targets need to be monitored to assess biodegradation. Microarray analysis may provide valuable insight into biodegradation of emerging contaminants, for which little is known regarding the microorganisms and degradation pathways involved. Microarrays have been used in research settings since 1996 but have only recently become commercially available for environmental applications. Microarrays have been used to document microbial diversity in a number of environments, including the petroleum release in the Gulf of Mexico (Hazen et al. 2010) and sites impacted by radionuclides like uranium (Chandler et al. 2010, Rastogi et al. 2010).

### What do microarrays do?

Environmental samples can contain thousands of different microorganisms and many different functional genes, some of which can serve as process-specific biomarkers. Phylogenetic microarrays evaluate community composition based on the presence/absence of microbial 16S rRNA genes present in a sample and answer the question, "Who is there?" A functional gene microarray targets genes involved in specific processes, for example a gene encoding a key enzyme involved in a degradation pathway, and can help answer questions about "potential activity." For example, functional gene microarray analysis can provide information on the capabilities of the microbial population to transform contaminants (e.g., degrade organic compounds, reduce metals such as Cr[VI]). Phylogenetic and functional microarrays can also be interrogated with RNA extracted from environmental samples and provide information about general activity (phylogenetic arrays) or about the activity of specific functional genes and pathways (functional array). Thus, microarrays can provide information about activity and determine, "Who is active?" and, "What pathway is active?"

### How are the data used?

The strength of the microarray approach is that many species or genes can be monitored simultaneously, and the overall responses of a microbial community to perturbations such as implementation of a remedy can be monitored over time or compared within impacted and background zones. A gram of soil or a liter of groundwater can contain billions of microorganisms, representing thousands of unique species that carry out different processes. Biodegradation of a contaminant of interest may require a single microbial population, a group of microorganisms, or a diverse community. In other words, a process of interest may be sufficiently monitored by looking at the dynamics of a couple of genes (e.g., genes encoding oxygenases involved in aerobic benzene biodegradation) in a couple of candidate species, while monitoring of more complex processes (e.g., nitrogen cycle, sulfur cycle, heavy-metals reduction) may be substantially improved by the analysis of hundreds or thousands of genes or assessment of flux of species present in a diverse community. Thus, microarrays can provide valuable insights for environmental remediation and monitoring.

### Example Environmental Remediation Questions Microarrays Can Help Answer

- **Site Characterization**
  - Assess current conditions and potential for biodegradation
    - What known microorganisms or functional genes are present?
    - How diverse is the microbial community?
  - Emerging contaminants
    - What types of microorganisms are present in impacted wells versus background wells?
- **Remediation**
  - Is monitored natural attenuation feasible?
    - Are supporting and contaminant-degrading microorganisms present in impacted zones?
    - Are competing microorganisms present?
  - Is biostimulation necessary? Should an amendment be added?
    - What kind of microorganisms will respond to the amendment?
- **Monitoring**
  - Monitored natural attenuation
    - What known microorganisms or functional genes are present?
    - Are competing microorganisms present?
    - How diverse is the microbial community?
    - Does microbial community structure (e.g., diversity) change over time?
    - What microorganisms are detected in impacted versus nonimpacted wells, and do these microorganisms support contaminant biodegradation?
  - Biostimulation
    - What changes were evident in the overall microbial community composition following amendment?
    - Did amendment addition promote growth of a specific group of microorganisms or increase of functional genes?
    - Is the shift in the microbial community consistent with the biostimulation strategy?
- **Closure**
  - How does the diversity of the microbial community in formerly impacted wells compare to the diversity of the microbial community in background wells?

### How do they work?

Microarrays are a collection of many short DNA strands, called “probes,” that are attached to a solid surface (e.g., a glass slide). The probes are selected for their specific, known DNA sequence, to which only complementary pieces of DNA (target) will bind (hybridize). After DNA is extracted from an environmental sample, it is fragmented and labeled with fluorescent chemicals and applied to the microarray. When hybridization (i.e., specific binding) occurs, the labeled DNA that complements its respective microarray probe is bound in place, producing characteristic fluorescent signals. DNA that does not have a complementary probe on the microarray slide is removed in a washing step. Detection and relative quantification are based on the fluorescent signal remaining after the washing step. This approach can also be applied to RNA obtained from the environmental sample. In this case, the RNA is transformed to complementary DNA (cDNA) in a step called “reverse transcription.” Hybridization of the cDNA to the array can provide information about activity. The strengths of the microarray approach are that many genes or species can be monitored simultaneously and the overall responses of a microbial community in response to remedial action can be monitored over temporal and spatial scales.

Figure 1 illustrates the results from a microarray analysis of DNA extracted from two individual environmental samples. For example, Sample A could be a groundwater sample obtained from a monitoring well in the contaminant source area, whereas Sample B could have been collected from a background well located upgradient of the site. The DNA from Sample A is labeled with a green fluorescent dye, and DNA from Sample B is labeled with a fluorescent red dye. Each position or “spot” in the microarray grid contains a specific gene probe. If Sample A contains complementary DNA target

sequences, the labeled DNA fragments will bind to the corresponding gene probes, producing a green signal at each of these positions. If Sample B DNA binds to the gene probes, the signals will be red. If DNA from both samples binds to the gene probes, a combination of both colors will appear (i.e., yellow). In the end, genes detected in Sample A only appear green, genes detected in Sample B only appear red, and genes detected in both samples appear yellow (a mix of green and red). Thus, in this example, the microarray results illustrate which microorganisms are unique to the contaminant source area (green), which are detected only in the background area (red), and which are present in both areas (yellow).




Figure 1. Example of results from a two-dye microarray.

### How are the data reported?

Phylogenetic microarray results are usually reported as a list of probable microorganisms (genus and species) detected in the sample. Similarly, functional gene microarray results include a list of the specific genes detected (e.g., a gene encoding nitrite reductase) and the gene type based on the biological process involved (e.g., denitrification). Statistical procedures have been developed that can aid in the interpretation of the results; however, microarray data interpretation requires expertise.

### Advantages

- Detection and relative quantification of thousands of organisms or functional genes in a single analysis.
- Information about gene expression (i.e., activity) can be obtained.
- Databases of known microorganisms and functional genes are becoming more comprehensive, making interpretation of results more meaningful and thus microarray analysis more applicable to environmental remediation. Microarrays provide a large quantity of information, which can be used to develop an understanding of the site that may not be possible using conventional environmental sampling and analytical testing. The microarray results may provide project managers with better information to use in the selection of remedial action alternatives or guide the selection of specific EMDs for efficient site monitoring. For example, microarrays can identify site-specific biomarker genes that provide meaningful information, and qPCR can then be applied to specifically monitor these genes.
- Gives an indication of the microbial diversity and possibly identifies the presence of microbes implicated in the biodegradation of the target contaminants. Microarrays can be based on both DNA and RNA, providing information on microbial community structure and metabolic activities, respectively.

## Limitations

- At the present time, few microarrays are commercially available that are relevant to environmental remediation.
- Careful design and thorough optimization and testing are needed to eliminate false positive signals (unspecific hybridization). Users should be sure to request documentation from the laboratory about the testing and validation of the microarrays.
- Quantification of the results can be difficult. Although recent studies have demonstrated relationships between signal intensity and target gene abundance, the dynamic range of the signal (i.e., the difference between the maximum and minimum signal) is limited and can hinder accurate quantification.
- Standardization of performance testing across different microarray platforms and guidelines for application and data interpretation are not readily available.
- The interpretation of data typically requires significant expertise, including knowledge of advanced statistical analyses.
- Microarray probes are based on genetic sequences of known microorganisms and biodegradation pathways cataloged in public databases—novel or as-yet undiscovered genes cannot be detected with microarrays. However, as new microorganisms and biodegradation pathways are identified, corresponding probes can be readily added to existing microarrays to expand the applicability of the technique to other contaminants and newly identified biodegradation pathways.

## Sampling Protocols

Almost any type of sample matrix (e.g., soil, sediment, groundwater) can be submitted for microarray analysis. Sampling usually involves collecting 10–20 g of soil or 1–2 L of groundwater and placement in sealed containers. Microarrays need a minimum of 2–5 µg of DNA; otherwise, it is necessary to amplify the sample prior to microarray analysis. The following items are typical requirements for microbiological sampling: (a) use of aseptic sample collection techniques and sterile containers, (b) shipment of the samples to the laboratory within 24 hours of sample collection, and (c) maintenance of the samples at 4°C during handling and transport to the laboratory. Sample collection techniques and containers may vary depending on the matrix sampled and the laboratory analyzing the samples. Users should work with the analytical laboratory to ensure sampling protocols for collecting, handling, and transporting the samples are in place and understood.

## Quality Assurance/Quality Control

To date, most EMDs do not have standardized protocols accepted by the U.S. Environmental Protection Agency (EPA) or other government agencies. However, EPA (2007) has an interim guidance for microarray analysis. In addition, most laboratories work under standard operating procedures (SOPs) and good laboratory practices, which can be provided to the user (e.g., consultant, state regulator) on request.

Currently, users can best ensure data quality by detailing the laboratory requirements in a site-specific quality assurance project plan (QAPP). This plan should include identification of the EMDs being used; the field sampling procedures, including preservation requirements; the SOPs of the laboratory performing the analysis; and any internal quality assurance/quality control information available (such as results for positive and negative controls). Specifically for microarrays, the arrays typically contain control probes and internal controls for analytical and technical performance of the system, as well as controls for normalization of signal. Standards currently exist for reporting data from microarray analysis (Brazma et al. 2001).

## Additional Information

He, Z., Y. Deng, J. D. Van Nostrand, Q. Tu, M. Xu, C. L. Hemme, X. Li, L. Wu, T. J. Gentry, Y. Yin, J. Liebich, T. C. Hazen, and J. Zho. 2010. "GeoChip 3.0 as a High-Throughput Tool for Analyzing

- Microbial Community Composition, Structure and Functional Activity," *The ISME Journal* **4**: 1167–79.
- He, Z., T. J. Gentry, C. W. Schadt, L. Wu, J. Liebich, S. C. Chong, Z. Huang, W. Wu, B. Gu, P. Jardine, C. Criddle, and J. Zhou. 2007. "GeoChip: A Comprehensive Microarray for Investigating Biogeochemical, Ecological and Environmental Processes," *The ISME Journal* **1**: 67–77. PMID 18043615. [www.ncbi.nlm.nih.gov/pubmed/18043615](http://www.ncbi.nlm.nih.gov/pubmed/18043615).
- U.S. Navy. 2009. "Molecular Biological Tools" website, vers. 1.1. [www.ert2.org/MBT/tool.aspx](http://www.ert2.org/MBT/tool.aspx).
- Wagner, M., H. Smidt, A. Loy, and J. Zhou. 2007. "Unraveling Microbial Communities with DNA-Microarrays: Challenges and Future Directions," *Microbial Ecology* **53**: 498–506. DOI 10.1007/s00248-006-9197-7.
- Wilson, K. H., W. J. Wilson, J. L. Radosevich, T. Z. DeSantis, V. S. Viswanathan, T. A. Kuczmarski, and G. L. Andersen, 2002. "High-Density Microarray of Small-Subunit Ribosomal DNA Probes," *Applied and Environmental Microbiology* **68**(5): 2535–41. PMID 11976131. [www.ncbi.nlm.nih.gov/pubmed/11976131](http://www.ncbi.nlm.nih.gov/pubmed/11976131).
- Wu, C. H., B. Sercu, L. C. Van De Werfhorst, J. Wong, T. Z. DeSantis, E. L. Brodie, T. C. Hazen, P. A. Holden, and G. L. Andersen. 2010. "Characterization of Coastal Urban Watershed Bacterial Communities Leads to Alternative Community-Based Indicators," *PLoS One* **5**(6): 11285. [www.plosone.org](http://www.plosone.org).
- Zhou, J. 2003. "Microarrays for Bacterial Detection and Community Analysis," *Current Opinion in Microbiology* **6**: 288–94.
- Zhou, J., Z. He, J. D. Van Nostrand, L. Wu, and Y. Deng. 2010. "Applying GeoChip Analysis to Disparate Microbial Communities," *Microbe* **5**(2): 60–65.

## References

- Biology Online. n.d. "Biology-Online Dictionary." [www.biology-online.org](http://www.biology-online.org).
- Brazma, A., P. Hingamp, J. Quackenbush, G. Sherlock, P. Spellman, C. Stoeckert, J. Aach, W. Ansorge, C. A. Ball, H. C. Causton, T. Gaasterland, P. Glenisson, F. C. P. Holstege, I. F. Kim, V. J. Makowitz, J. C. Matese, H. Parkinson, A. Robinson, U. Sarkans, S. Schulze-Kremer, J. Stewart, R. Taylor, J. Vilo, and M. Vingron. 2001. "Minimum Information about a Microarray Experiment (MIAME): Toward Standards for Microarray Data," *Nature Genetics* **29**: 365–71.
- Chandler, D. P., A. Kukhtin, R. Mokhiber, C. Knickerbocker, D. Ogles, G. Rudy, J. Golova, P. Long, and A. Peacock. 2010. "Monitoring Microbial Community Structure and Dynamics during In Situ U(VI) Bioremediation with a Field-Portable Microarray Analysis System," *Environmental Science and Technology* **44**: 5516–22. PMID 20560650. [www.ncbi.nlm.nih.gov/pubmed/20560650](http://www.ncbi.nlm.nih.gov/pubmed/20560650).
- EPA (U.S. Environmental Protection Agency). 2007. *Interim Guidance for Microarray-Based Assays*. External review draft. [www.epa.gov/spc/pdfs/epa\\_interim\\_guidance\\_for\\_microarray-based\\_assays-external-review\\_draft](http://www.epa.gov/spc/pdfs/epa_interim_guidance_for_microarray-based_assays-external-review_draft).
- EPA. 2011. "Glossary of Technical Terms." [www.epa.gov/oust/cat/tumgloss.htm#a](http://www.epa.gov/oust/cat/tumgloss.htm#a).
- Hazen, T. C., E. A. Dubinsky, T. Z. DeSantis, G. L. Andersen, Y. M. Piceno, N. Singh, J. K. Jansson, A. Probst, S. E. Borglin, J. L. Fortney, W. T. Stringfellow, M. Bill, M. E. Conrad, L. M. Tom, K. L. Chavarria, T. R. Alusi, R. Lamendella, D. C. Joyner, C. Spier, J. Baelum, M. Auer, M. L. Zemla, R. Chakraborty, E. L. Sonnenthal, P. D'haeseleer, H. N. Holman, S. Osman, Z. Lu, J. D. Van Nostrand, Y. Deng, J. Zhou, and O. U. Mason. 2010. "Deep-Seas Oil Plume Enriches Indigenous Oil-Degrading Bacteria," *Science* **330**(6001): 204–08.
- Nannipieri, P., J. Ascher, M. T. Ceccherini, L. Landi, G. Pietramellara, and G. Renella. 2003. "Microbial Diversity and Soil Functions," *European Journal of Soil Science* **54**: 655–70.
- Rastogi G, S. Osman, P. A. Vaishampayan, G. L. Andersen, L. D. Stetler, and R. K. Sani. 2010. "Microbial Diversity in Uranium Mining-Impacted Soils as Revealed by High-Density 16S Microarray and Clone Library," *Microbial Ecology* **59**(1): 94–108. PMID 19888627. [www.ncbi.nlm.nih.gov/pubmed/19888627](http://www.ncbi.nlm.nih.gov/pubmed/19888627).



## Glossary

**16S rRNA**—A subunit of the ribosome composed of ribonucleic acid (RNA). The RNA sequence is used to classify and identify microorganisms (e.g., genus and species).

**biodegradation**—A process by which microorganisms transform or alter (through metabolic or enzymatic action) the structure of chemicals introduced into the environment (EPA 2011).

**biomarker**—A distinctive (unique) characteristic of a biomolecule that can be measured and used as an indicator of a target microorganism or biological process. For example, a specific DNA sequence (used as a probe on a microarray) could be a biomarker for a particular microorganism (e.g., *Desulfotomaculum*).

**functional gene**—A segment of DNA that encodes an enzyme or other protein that performs a known biochemical reaction. For example, the functional gene *tceA* encodes the reductive dehalogenase enzyme that initiates reductive dechlorination of trichloroethene. Other genes can code for RNA entities which can regulate the activity of other DNA target sequences.

**genus**—A category of organism classification (taxonomy). A particular genus is a group of related species. For example, *Pseudomonas* is a genus of bacteria.

**microarray probe**—A short, defined segment of DNA designed to bind with the target gene if found in the environmental sample. The probe is attached to the solid surface of the microarray.

**microbial community**—The microorganisms present in a particular sample.

**microbial diversity**—Microbial diversity can have many definitions but in this context generally refers to the number of different microbial species and their relative abundance in an environmental sample (Nannipieri et al. 2003).

**nitrite reductase gene**—Functional genes encoding the enzymes that catalyze nitrite reduction. Nitrite reductase genes are commonly used as the target gene to detect microorganisms capable of denitrification.

**phylogeny (phylogenetic analysis)**—Classification of microorganisms into groups (e.g., genus and species) based in part on the rRNA sequences.

**ribosome**—A multicomponent biological molecule which is part of the protein-synthesizing machinery of the cell.

**species**—The lowest taxonomic rank and the most basic unit or category of biological classification (Biology Online n.d.).

### EMD Team Contact

Robert Mueller, Team Leader

New Jersey Department of Environmental Protection

[bob.mueller@dep.state.nj.us](mailto:bob.mueller@dep.state.nj.us), (609) 984-3910



ITRC is affiliated with the  
Environmental Council of the States

### ACKNOWLEDGEMENTS

The members of the Interstate Technology & Regulatory Council (ITRC) Environmental Molecular Diagnostics (EMD) Team wish to acknowledge the individuals, organizations, and agencies that contributed to this set of fact sheets.

As part of the broader ITRC effort, the EMD Team effort is funded by the U.S. Department of Energy, U.S. Department of Defense, and the U.S. Environmental Protection Agency and through ITRC's Industry Affiliates Program.

The EMD Team wishes to thank the ITRC external reviewers and the peer reviewers who contributed comments and suggestions that were of great help to the team in finalizing the fact sheets. The EMD Team also wishes to recognize and thank Bonnie Pierce, formerly of the Wyoming Department of Environmental Quality, who was our team co-leader during 2010 and whose leadership helped guide the development of these fact sheets.

The EMD Team worked hard to develop, review, and revise this set of fact sheets. The team recognizes the great value of teamwork and thanks everyone who participated—named and unnamed, ITRC staff, ITRC Point of Contact, or team member.

The EMD Team recognizes the efforts and important contributions of the following state environmental personnel: James Fish, Alaska Department of Environmental Conservation; Christine Brown, Vivek Mathrani, Sara Michael, and Claudio Sorrentino, California Department of Toxic Substance Control; Cleet Carlton, California Regional Water Quality Control Board; Leslie Smith, Florida Department of Environmental Protection; Amanda Howell and Undine Johnson, Georgia Environmental Protection Division; Robert Mueller, New Jersey Department of Environmental Protection, EMD Team Leader; and Ramesh Belani, Pennsylvania Department of Environmental Protection.

The EMD Team recognizes the efforts and valuable contributions of the following stakeholder and academic representatives: Peter Strauss, PM Strauss & Associates; Michael Hyman, North Carolina State University; Frank Löffler, University of Tennessee; Paul Philp, University of Oklahoma; Kerry Sublette, University of Tulsa; and Jennifer Weidhaas, West Virginia University.

The EMD Team recognizes the efforts and valuable contributions of the following federal personnel: Adria Bodour and John Gillette, AFCEE; Ann Miracle, DOE, Pacific Northwest National Laboratory; Hans Stroo, SERDP/ESTCP; Cheryl A. Hawkins and Ann Keeley, U.S. EPA; and Carmen Lebrón, U.S. Navy.

The EMD Team recognizes the efforts and valuable contributions of the following consultants and industry representatives: Stephen Koenigsberg, Adventus Americas, Inc.; Rebecca Mora, Chad Roper, Matthew Mesarch, and Jing Zhou, AECOM Environment; Jessica Goin, Anchor QEA; Caitlin Bell, Rula Deeb, and Denise Nelson, ARCADIS; Ramona Darlington, Battelle Memorial Institute; Stephanie Fiorenza, BP; M. Hope Lee, Tamzen Macbeth, and Ryan Wymore, CDM; David Duncklee, Duncklee and Dunham; William Berti, DuPont; Eric Raes, Engineering and Land Planning Associates, Inc.; Devon Rowe, ENVIRON; David Major and Erik Petrovskis, Geosyntec Consultants; Ioana Petrisor, Haley & Aldrich, Inc.; Sophia Drugan, Kleinfelder, Inc.; Brett Baldwin, Dora Ogles, and Greg Davis, Microbial Insights, Inc.; Pat McLoughlin Microseeps, Inc.; Lesley Hay Wilson, Sage Risk Solutions, LLC; and Paul Hatzinger, Shaw Environmental.

#### **ABOUT ITRC**

The Interstate Technology & Regulatory Council (ITRC) is a public-private coalition working to reduce barriers to the use of innovative environmental technologies and approaches so that compliance costs are reduced and cleanup efficacy is maximized. ITRC produces documents and training that broaden and deepen technical knowledge and expedite quality regulatory decision making while protecting human health and the environment. With private- and public-sector members from all 50 states and the District of Columbia, ITRC truly provides a national perspective. More information on ITRC is available at [www.itrcweb.org](http://www.itrcweb.org).

ITRC is a program of the Environmental Research Institute of the States (ERIS), a 501(c)(3) organization incorporated in the District of Columbia and managed by the Environmental Council of the States (ECOS). ECOS is the national, nonprofit, nonpartisan association representing the state and territorial environmental commissioners. Its mission is to serve as a champion for states; to provide a clearinghouse of information for state environmental commissioners; to promote coordination in environmental management; and to articulate state positions on environmental issues to Congress, federal agencies, and the public.

#### **DISCLAIMER**

This material was prepared as an account of work sponsored by an agency of the U.S. Government. Neither the U.S. Government nor any agency thereof, nor any of their employees, makes any warranty, express or implied, or assumes any legal liability or responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise does not necessarily constitute or imply its endorsement, recommendation, or favoring by the U.S. Government or any agency thereof. The views and opinions of authors expressed herein do not necessarily state or reflect those of the U.S. Government or any agency thereof, and no official endorsement should be inferred.

The information provided in documents, training curricula, and other print or electronic materials created by the Interstate Technology & Council ("ITRC Products") is intended as a general reference to help regulators and others develop a consistent approach to their evaluation, regulatory approval, and deployment of environmental technologies. The information in ITRC Products is formulated to be reliable and accurate. However, the information is provided "as is," and use of this information is at the users' own risk.

ITRC Products do not necessarily address all applicable health and safety risks and precautions with respect to particular materials, conditions, or procedures in specific applications of any technology. Consequently, ITRC recommends consulting applicable standards, laws, regulations, suppliers of materials, and material safety data sheets for information concerning safety and health risks and precautions and compliance with then-applicable laws and regulations. ITRC, ERIS, and ECOS shall not be liable in the event of any conflict between information in ITRC Products and such laws, regulations, and/or other ordinances. ITRC Product content may be revised or withdrawn at any time without prior notice.

ITRC, ERIS, and ECOS make no representations or warranties, express or implied, with respect to information in ITRC Products and specifically disclaim all warranties to the fullest extent permitted by law (including, but not limited to, merchantability or fitness for a particular purpose). ITRC, ERIS, and ECOS will not accept liability for damages of any kind that result from acting upon or using this information.

ITRC, ERIS, and ECOS do not endorse or recommend the use of specific technologies or technology providers through ITRC Products. Reference to technologies, products, or services offered by other parties does not constitute a guarantee by ITRC, ERIS, and ECOS of the quality or value of those technologies, products, or services. Information in ITRC Products is for general reference only; it should not be construed as definitive guidance for any specific site and is not a substitute for consultation with qualified professional advisors.



## Regulatory Acceptance for New Solutions

Documents, free Internet-based training, contact information

[www.itrcweb.org](http://www.itrcweb.org)